Recent advances in picosecond and hybrid femtosecond/picosecond CARS techniques for gas phase analysis

Christopher J. Kliewer, Combustion Research Facility, Sandia National Laboratories, CA USA

Accurate time- and spatially resolved measurement of temperature remains a critical focus of combustion diagnostics. The local temperature field and gradient not only govern chemical reaction rates but also physical quantities such as gas expansion and heat transfer. Coherent anti-Stokes Raman spectroscopy (CARS) has been applied for concentration measurements and ro-vibrational thermometry in gas-phase applications for more than three decades, and is often held as the gold standard optical technique for the non-intrusive determination of temperature in gas-phase systems. Current developments in picosecond (ps) coherent anti-Stokes Raman spectroscopy (CARS) and hybrid femtosecond/picosecond (fs/ps) CARS techniques have enabled breakthroughs in measurement capabilities for gas phase molecular analysis.

Utilizing a ps probe pulse in gas phase CARS investigations is in many ways an optimized approach. Collisional dephasing of many gas phase molecules occurs on the ps timescale for a range of pressures, thus, ps probe pulses are uniquely formatted to directly probe such processes. For instance, time-resolved ps-CARS has been recently employed to acquire high-fidelity measurements of collisional dephasing and, hence, Raman broadening coefficients [1-4], although for quickly dephasing molecules care must be taken to quantify the time-convolution with the probe window. A key advantage of ps-probing in CARS experiments is the capability for both excellent frequency resolution (transform-limited bandwidth a small fraction of 1 cm⁻¹), while simultaneously maintaining the advantages of time-resolved CARS with high peak power pulses, such as very high signal level, suppression of nonresonant background, and direct measurement of coherence dephasing effects.

Hybrid fs/ps CARS techniques [5] take advantage of *both* impulsive fs coherence excitation and high power, well-resolved, ps probing. Many breakthrough developments have recently been enabled by further developments with hybrid fs/ps CARS. Two-beam 1D ultrabroadband CARS has been demonstrated for, in principle, the simultaneous measurement of all Raman-active species and molecular manifolds in a gas phase mixture [6], demonstrated in Fig. 1. The first gas phase 2D-CARS measurements generating thousands of spatially correlated CARS spectra within a single laser shot have been demonstrated [7] in the harsh environment of a flame. A very recent development within our laboratory has been the development of ultrabroadband 2D-CARS – enabling planar imaging of Raman active-modes and spectra with an excitation bandwidth demonstrated beyond 4000 cm⁻¹.




Fig. 1. (a) CCD image of a single shot acquisition of ultrabroadband CARS in a mixture of N_2 , H_2 , and methane. The log of signal intensity is presented. (b) CARS spectrum taken from a single pixel-row of the CCD image in panel (a). Error bars represent single-shot standard deviation in signal intensity when normalized to the N_2 pure-rotational branch.

1. W. D. Kulatilaka, P. S. Hsu, H. U. Stauffer, J. R. Gord, and S. Roy, Appl. Phys. Lett. 97 (2010).

2. C. J. Kliewer, Y. Gao, T. Seeger, J. Kiefer, B. D. Patterson, and T. B. Settersten, Proc. Combust. Inst. **33**, 831-838 (2011).

3. C. Kliewer, A. Bohlin, E. Nordström, B. Patterson, P.-E. Bengtsson, and T. Settersten, Appl. Phys. B. **108**, 419-426 (2012).

4. S. Roy, P. S. Hsu, N. Jiang, J. R. Gord, W. D. Kulatilaka, H. U. Stauffer, and J. R. Gord, J. Chem. Phys. **138**, 024201 (2013).

5. J. D. Miller, M. N. Slipchenko, T. R. Meyer, H. U. Stauffer, and J. R. Gord, Opt. Lett. **35**, 2430-2432 (2010).

- 6. A. Bohlin, and C. J. Kliewer, Appl. Phys. Lett. **104**, 031107 (2014).
- 7. A. Bohlin, and C. J. Kliewer, The Journal of Physical Chemistry Letters, 1243-1248 (2014).